

**Abutment-free bone-anchored hearing devices in children: Initial results and experience**

Published *Int J Pediatr Otorhinolaryngol.* 2014 May

Summary of a Clinical Study

Sri Kiran Chennupati,

principal investigator

<b>Clinic/ Physicians</b>	<ul style="list-style-type: none"> <li>Section of Pediatric Otolaryngology, Drexel University College of Medicine, St. Christopher's Hospital for Children, 3601 A Street, Philadelphia, PA 19134, USA</li> </ul>														
<b>Objective</b>	<ul style="list-style-type: none"> <li>Aaron Centric, <b>Sri Kiran Chennupati</b></li> </ul>														
<b>Patients and Study Design</b>	<ul style="list-style-type: none"> <li>Assessing the Sophono Alpha 1, a new transcutaneous abutment-free system, as compared with other percutaneous bone conduction systems, which have a complication rate of more than 20% that arise mainly from the skin-penetrating abutment.</li> </ul>														
<b>Evaluations</b>	<ul style="list-style-type: none"> <li>5 patients – all underwent Sophono implantation</li> </ul>	<ul style="list-style-type: none"> <li>The patients were between 5 years to 12 years, median 10 years</li> <li>1 patient moderate bilateral conductive hearing loss, 1 patient unilateral sensorineural hearing loss, 2 patients had moderate conductive hearing loss, 1 patient severe conductive hearing loss</li> <li>Average length of follow up: 145 days</li> </ul>	<ul style="list-style-type: none"> <li>Conduction of a retrospective chart review of the first five patients who underwent implantation with the Sophono abutment-free bone conduction hearing system with the Alpha 1 processor at the institution and reporting on the patients' pre- and post-operative audiometric data and clinical courses.</li> </ul>												
<b>Other Parameters</b>	<table border="1"> <tr> <td>1 week after implantation</td> <td>Wound check</td> </tr> <tr> <td>6 weeks after implantation</td> <td>Individual frequency data using conventional audiometry, pure-tone averages (PTAs), speech recognition thresholds (SRT), and masking of the unaided ear.</td> </tr> </table>	1 week after implantation	Wound check	6 weeks after implantation	Individual frequency data using conventional audiometry, pure-tone averages (PTAs), speech recognition thresholds (SRT), and masking of the unaided ear.	<table border="1"> <tr> <td>Average length of time from implantation to activation</td> <td><b>145 days (SD 57)</b></td> </tr> <tr> <td>Sophono Processor Type</td> <td><b>Alpha 1</b></td> </tr> </table>	Average length of time from implantation to activation	<b>145 days (SD 57)</b>	Sophono Processor Type	<b>Alpha 1</b>					
1 week after implantation	Wound check														
6 weeks after implantation	Individual frequency data using conventional audiometry, pure-tone averages (PTAs), speech recognition thresholds (SRT), and masking of the unaided ear.														
Average length of time from implantation to activation	<b>145 days (SD 57)</b>														
Sophono Processor Type	<b>Alpha 1</b>														
<b>Timeline</b>	<ul style="list-style-type: none"> <li>At the time of submission the average length of follow-up was 145 days (SD 57)</li> </ul>														
<b>Clinical Results</b>	<table border="1"> <thead> <tr> <th></th> <th>Pre-Op (unaided)</th> <th>Post-Op (aided)</th> <th>Change (improvement)</th> </tr> </thead> <tbody> <tr> <td>Pure-Tone Average</td> <td>57 dB</td> <td>25 dB</td> <td>32 dB</td> </tr> <tr> <td>SRT</td> <td>56 dB</td> <td>28 dB</td> <td>28 dB</td> </tr> </tbody> </table>				Pre-Op (unaided)	Post-Op (aided)	Change (improvement)	Pure-Tone Average	57 dB	25 dB	32 dB	SRT	56 dB	28 dB	28 dB
	Pre-Op (unaided)	Post-Op (aided)	Change (improvement)												
Pure-Tone Average	57 dB	25 dB	32 dB												
SRT	56 dB	28 dB	28 dB												
<b>Conclusion</b>	<ul style="list-style-type: none"> <li>Pressure points over the implant for one patient were resolved after compliance to the proper post-activation instructions as provided by Sophono, by adjusting the appropriate strength of the magnet and short spans of time of use until full accommodation.</li> <li>All patients were responding in the normal to mild hearing loss range in the operated ear after device activation. Average improvement across individual frequencies was between 17 and 37 dB (SD 5.5–11 dB).</li> <li>The surgical procedure is relatively simple, and morbidity is low while recovery is fast. Audiometric results have been consistent with published data on other bone-anchored hearing devices.</li> </ul>														

Rx only. Refer to product instruction manual/package insert for instructions, warnings, precautions and contraindications.

For further information, please call Medtronic ENT at 800.874.5797 or consult Medtronic's website at [www.medtronicent.com](http://www.medtronicent.com).



Now owned by

## Medtronic

### Sophono Inc.

5744 Central Avenue  
Suite 100  
Boulder, CO 80301

### International Telephone Numbers

Adriatic Region 385-1-488-1120	Israel 972-9-972-4400
Argentina 54-11-4898-5700	Italy 39-02-24137-324
Australia 1-800-668-670	Japan 81-3-6430-2017
Baltic Region 37-1-67560226	Korea 82-2-3404-3600
Belgium 32-2456-09-09	Lebanon 961-1-370-670
Brazil 55-11-2182-9200	Luxembourg 32-2456-09-09
Canada 1-800-217-1617	Malaysia 60-37-953-4800
Chile 56-2-2655-5110	Mexico 52-55-11-02-90-30
China 86-21-20325888	Netherlands 31-45-566-8800
Colombia 57-1-742-7300	Poland 48-22-4656900
Czech Republic 420-2-9657-9580	Russian Fed. 7-495-580-73-77
France 33-470-679-800	Singapore 65-6776-6255
Germany 49-2159-8149-209	South Africa 27-11-466-1820
Greece 30-210-67-79-099	Spain 34-91-625-05-40
Hong Kong 852-2919-1300	Taiwan 886-2-2183-6000
Hungary 36-188-90600	Thailand 662-232-7400
India 91-22-26836733	UK 44-1923-205-166

Sophono™ and *The Leader in Magnetic Bone Conduction Hearing* are registered trademarks of Medtronic, Inc.  
© 2016, Medtronic, Inc. S0584-08 Rev B UC201700050 EN 07.2016